



VERTICAL MINI GEL SYSTEM

BIO 05-01

Complete Solution for Life Science Research

VERTICAL MINI GEL SYSTEM (BIO 05-01)



Technical Specifications

Unit Dimensions (L x W x H)	15 x 18 x 15cm
Gel Dimensions (L x W)	8 x 7 cm (1 Gel)
Buffer Volume	100 + 150 ml
Maximum Sample Capacity	7
Combs	7 Well Teflon Comb 0.5mm,1mm,1.5mm – 1 No Each
Spacers	Teflon Spacers 0.5mm,1mm,1.5mm – 2 No's Each
Platinum Electrodes	Red and Black (1 Each)
Glass Plates	Notched and Rectangular (2 Sets)
Gasket	Fixed
Clamps and Screws	1 Set
Leveling Screws	3 No's
Connecting Cord	1 Pair
Sprit Level(Fixed)	1 No
User Manual	

BIO WORLD

1247, 6th Main, 5th Stage, 1st Phase, B.E.M.L Layout , R.R Nagar, Bangalore – 560098

Website: www.bioworldproducts.com ; Email : info@bioworldindia.in

Mobile: +91-9591779325,7760658542,9448351365