



VERTICAL MAXI GEL SYSTEM

BIO 05-05

Complete Solution for Life Science Research

VERTICAL MAXI GEL SYSTEM (BIO 05-05)



Technical Specifications

Unit Dimensions (L x W x H)	15 x 18 x 15cm
Gel Dimensions (L x W)	16 x 20 cm (1 Gel)
Buffer Volume	150 + 200 ml
Maximum Sample Capacity	20
Combs	20 Well Teflon Comb 0.5mm,1mm,1.5mm – 1 No Each
Spacers	Teflon Spacers 0.5mm,1mm,1.5mm – 2 No's Each
Platinum Electrodes	Red and Black (1 Each)
Glass Plates	Notched and Rectangular (2 Sets)
Gasket	Fixed
Clamps and Screws	1 Set
Leveling Screws	3 No's
Connecting Cord	1 Pair
Sprit Level (Fixed)	1 No

User Manual

BIO WORLD

1247, 6th Main, 5th Stage, 1st Phase, B.E.M.L Layout , R.R Nagar, Bangalore – 560098

Website: www.bioworldproducts.com ; Email : info@bioworldindia.in

Mobile: +91-9591779325,7760658542,9448351365