



DIGITAL HOT PLATE

Complete Solution for Life Science Research



Technical Specification:

Description	Digital Hot Plate
Top Plate	Made of thick plate SS 304grade
Outer Chamber	Outer chamber made of mild steel and finished with powder coated.
Insulation	Top Plate insulated with high grade Ceramic Fibre
Heaters	Coiled Air Heaters are made of high grade Kanthal wires of suitable wattage.
Temperature Range	5 °C above ambient to 300°C
Temperature Controller / Resolution	Microprocessor based PID temperature controller with dual display (Process value and set value), low and high limit setting with Cr Al thermocouple and Timer control. / 0.1°C.
Power Supply	230 Volts Single Phase 50Hz, AC supply.
User Manual	Yes

Ordering Information

Cat#	Description	Qty
BIO HP-01	Hot Plate 10" x 10"	1 No
BIO HP-02	Hot Plate 12" x 12"	1 No
BIO HP-03	Hot Plate 12" x 18"	1 No
BIO HP-04	Hot Plate 18" x 24"	1 No

BIO WORLD

1247, 6th Main, 5th Stage, 1st Phase, B.E.M.L Layout , R.R Nagar, Bangalore – 560098

Website: www.bioworldproducts.com ; Email : info@bioworldindia.in

Mobile: +91-9591779325,7760658542,9448351365