



GEL ROCKER

BIO GR-01

Complete Solution for Life Science Research

Gel Rocker (#BIO GR-01)



Gel Rocker can be used for many applications like gel staining, culture aeration, bathing, blotting membranes etc. A metal platform with non-slip rubber mat moves in a gentle 'See Saw' Motion. The Movement is controlled precisely with a DC motor. Cushion and rubber feet prevent unit from Slipping. Individual platform rocker 2 Samples can be stained at a time.



Technical Specification

Description	Gel Rocker (#BIO GR-01)
Cabinet	MS Powder Coated with SS Working Top
Dimensions	30 x 30 x 11.5 cm
Platform Size	25 x 25 cm
Weight	4 Kg
Movement	See Saw
Operating Voltage	230v ac@ 50 Hrtz
User Manual	Yes

Ordering Information

Cat#	Description	Qty
BIO-GR-01	Gel Rocker (Manual) – Single Platform	1 No
BIO-GR-02	Gel Rocker (Timer) – Single Platform	1 No
BIO-GR-03	Gel Rocker (Timer) – Dual Platform	1 No

BIO WORLD

1247, 6th Main, 5th Stage, 1st Phase, B.E.M.L Layout , R.R Nagar, Bangalore – 560098

Website: www.bioworldproducts.com ; Email : info@bioworldindia.in

Mobile: +91-9591779325,7760658542,9448351365