



INCUBATOR ORBITAL SHAKER FLOOR MODEL IOS-03

Complete Solution for Life Science Research

Incubator Orbital Shaker (Floor Model)



Incubator Orbital shaker is in almost every kind of laboratory. It is a very important device in order to perform tests and researches. The use of this instrument is very common; you can see them on a workspace, on top of lab's bench tops growing cultures. Selecting the right shaker is a very important step when it comes to building your laboratory.

Mixing samples, solutions and different substances is a very common technique in most medical, biological and chemistry labs. That is why we are presenting you with this orbital shaker model. Being this a mandatory instrument in your workspace.

Technical Specification

Description	Incubator Orbital Shaker Floor Model
Outer Chamber	Aluminum powder coated
Inner Chamber	Aluminum powder coated
Speed	30-250 rpm
Motion:	Orbital
Control:	Feedback control System
Temperature range:	Above ambient to 60°C
Controller:	Digital temperature controller PID type
Accuracy:	+/- 0.5°C
Flask Holder	UNIVERSAL TRAY WITH SPRINGS
Overall Dimension:	W x D x H (730 x 750 x 770)mm
Capacity:	500ml x 16nos

Ordering Information

Cat#	Description	Qty
BIO-IOS-04	Incubator Orbital Shaker Cooling 450 x 450 mm	1 No

BIO WORLD

1247, 6th Main, 5th Stage, 1st Phase, B.E.M.L Layout , R.R Nagar, Bangalore – 560098

Website: www.bioworldproducts.com ; Email : info@bioworldindia.in

Mobile: +91-9591779325,7760658542,9448351365