



ORBITAL SHAKER PROGRAMMABLE Microprocessor Control – LCD Display

Complete Solution for Life Science Research

Orbital Shaker Programmable (#BIO OS-02)



Orbital shaker in almost every kind of laboratory. It is a very important device in order to perform tests and researches. The use of this instrument is very common; you can see them on a workspace, on top of lab's bench tops growing cultures. Selecting the right shaker is a very important step when it comes to building your laboratory.

Mixing samples, solutions and different substances is a very common technique in most medical, biological and chemistry labs. That is why we are presenting you with this orbital shaker model. Being this a mandatory instrument in your workspace.

Technical Specification

Description	Orbital Shaker Programmable (# BIO OS-02)
Cabinet	GI Body with Powder Coated
Speed	40 – 250rpm
Display	Jumbo LCD Display
Control Type	PID - RPM Controller
Platform Holder	Spring Type (Universal)
Movement	Orbital Motion
Speed	40 to 250 Rpm
Timer	Yes
Program	Can set up to 20 programs
Operating Voltage	AC 220V, 50 Hz
User Manual	Yes

Ordering Information

Cat#	Description	Qty
BIO-OS-01	Orbital Shaker (Electronic)	1 No
BIO-OS-02	Orbital Shaker (Electronic – Programmable)	1 No

BIO WORLD

1247, 6th Main, 5th Stage, 1st Phase, B.E.M.L Layout , R.R Nagar, Bangalore – 560098

Website: www.bioworldproducts.com ; Email : info@bioworldindia.in

Mobile: +91-9591779325,7760658542,9448351365