



GEL DOCUMENTATION
COLLEGE MODEL
BIO-GD-01

Complete Solution for Life Science Research

GEL DOCUMENTATION COLLEGE MODEL (#BIO-GD-01)



Technical Specification

Description	Gel Documentation College (#BIO-GD-01)
Material	MS with Powder Coated
Camera	Digital DSLR Camera 12 MP
Power Cords	For UV and Camera Separately
Camera Features	Capture the gel and viewing
Camera Hood	Doom Type
UV Transilluminator	Separately Provided
Viewing Area	21 x 21 Cms
Wavelength	312nm
Detection Limit	0.1ng
Application	DNA and RNA gels imaging
User Manual	Yes
Assembly	UV transilluminator , DARK HOOD and camera needs to be fixed separately

Ordering Information:

Cat#	Description	Qty
BIO-GD-01	Gel Documentation (College Model)	1 No
BIO-GD-02	Gel Documentation (Research) CMOS with 1D Software	1 No
BIO-GD-03	Gel Documentation (Research) CCD with 1D Image Analysis Software	1 No

BIO WORLD

1247, 6th Main, 5th Stage, 1st Phase, B.E.M.L Layout , R.R Nagar, Bangalore – 560098

Website: www.bioworldproducts.com ; Email: info@bioworldindia.in

Mobile: +91-9591779325,7760658542,9448351365